Lichthärtbares ONE-FOR-MOST SHADES Composite Restaurationsmaterial

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Light Cure ONE-FOR-MOST SHADES composite restorative

LC CHROME is a light cure radiopaque composite for anterior and posterior restorations. LC CHROME can be used for all natural tooth shades due to its continuous color matching properties. Beside its permanent gloss and toothlike opalescence LC CHROME shows excellent physical properties and good handling features. High gloss polishability demands only little effort. LC CHROME contains dimethacrylate based resins and inorganic fillers < 0.2 µm. The total filler content is 80 % by weight and 75 % by volume. Radiopacity of LC CHROME is 200 % aluminium (radiopacity of enamel: ~200 % aluminium, of dentine: ~100 % aluminium).

LC CHROME meets the requirements of **ISO 4049**, type 1, class 2, group 1.

Indications LC CHROME

- Restorations of all cavity classes
- Direct veneers
- Diastema closure
- Repair of ceramic restorations (e.g. veneers, crowns)

Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with pulp capping material (e.g. LC ResinCal PC or CALCIUM HYDROXIDE LINER).

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Preparation of Tooth

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water).

Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

For pulp protection areas close to the pulp should be covered with a with pulp capping material (e.g. LC ResinCal PC or CALCIUM HYDROXIDE LINER). Refer to the corresponding user instructions.

Apply a common adhesive for enamel and dentine bonding according to the corresponding manufacturer instructions.

The application of the self-etching adhesive LC ETCHBOND X1 is recommended. Follow the instructions for LC ETCHBOND X1.

Application and Curing of LC CHROME

In case of deep and extended cavities best results are obtained with an application of a base liner using a light cure resin modified cement or a light cure flowable composite, followed by the application of LC CHROME.

In case of small cavities LC CHROME can be applied directly.

Note:

In case of cavities of dark or discolored teeth use a light cure blocker to ensure that the restoration may not appear dark.

Apply LC CHROME in **layers of max. 2 mm**. Light cure each increment separately. LC CHROME is easy to model. Place it by using a suitable instrument. Matrix strips may be used.

By using a polymerization unit (wavelength range 400-500 nm) with a **light intensity of at least 1000mW/cm²**, cure each increment separately for **20 seconds**.

Hold the light emission window as close as possible to the filling material.

Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burs.

Polish to high gloss with polishing discs or polishing brushes.

Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

Additional Notes/Warnings

- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

Warranty

First Scientific Dental GmbH warrants this product will be free from defects in material and manufacture. First Scientific Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and First Scientific Dental GmbH's sole obligation shall be repair or replacement of the First Scientific Dental GmbH product.

Limitation of Liability

Except where prohibited by law, First Scientific Dental GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!

For dental use only!