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Self adhesive light cure flowable composite

LC FLOWFILL SE is a **self adhesive light cure flowable composite for anterior and limited posterior restorations**. Beside its good mechanical properties, the physical properties of LC FLOWFILL SE are comparable to conventional flowable composites. The optimized flow behavior enables easy application of the material.

The self adhesives properties of LC FLOWFILL SE require no etching, priming and bonding of enamel and dentine in the restorative therapy. For application as pit and fissure sealant, enamel etching for cleaning the surface is recommended.

LC FLOWFILL SE is based on methacrylates and inorganic fillers and inorganic filler particles of 0.05-1 µm. The total filler load is 59% and the total filler volume 45%. Radiopacity of LC FLOWFILL SE is 150 % aluminium (radiopacity of enamel: ~200 % aluminium, of dentine:~ 100 % aluminium). The composite material meets the requirements of **DIN EN ISO 4049**, type 1, class 2, group 1.

Indications

- Small class I restorations
- Lining of class I and II cavities
- Pit and fissure sealing
- Extended fissure sealing
- Filling of undercuts

Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol can impair the hardening of the composite and causes discoloration. Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow to come into contact with materials containing eugenol.

Preparation of Tooth and Color Matching

Clean the tooth with a fluoride-free polishing paste prior to preparation and color matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate LC FLOWFILL SE shade.

LC FLOWFILL SE is shaded according to VITA®-Shades.

Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction. If possible, margins should have a slight **(0.5 - 1.0 mm)** bevel placed in the enamel to increase the surface area for greater bond strength.

For pulp protection areas close to the pulp should be covered with a thin layer of e.g. LC ResinCal PC or CALCIUM HYDROXIDE LINER. Please refer to the instructions of the product used.

Placement of LC FLOWFILL SE

Remove the cap from the LC FLOWFILL SE-syringe and attach a NeedleTip. Squeeze out the first small drop on a pad to ensure safe handling of LC FLOWFILL SE.

After use discard the NeedleTip and attach the cap to the syringe again. NeedleTips are for single use only due to hygienic reasons.

Small class I restorations and lining of class I and II

Place LC FLOWFILL SE with the NeedleTip in a thin layer directly onto the prepared areas. With the supplied microbrush work it into the entire cavity walls and beveled areas under slight pressure for 20s. Remove excess material around the margins with a microbrush.

Cure LC FLOWFILL SE as follows:

Lighter shades (e.g. A1, A2, B1) **20 seconds**

Darker shades (e.g. A3) **30 seconds**

After lining of the cavity walls and the beveled areas build the restoration with LC FLOWFILL SE in **layers of max. 2 mm**.

Light cure each increment separately.

Lighter shades (e.g. A1, A2, B1) **20 seconds**

Darker shades (e.g. A3) **30 seconds**

Pit and fissure sealing

Clean the enamel with a fluoride-free paste and rinse thoroughly. Dry in a water and oil-free air stream leaving a slightly moist shining surface on the enamel.

Note:

In addition to cleaning as described above, etching of enamel with a 37 % phosphoric acid etching gel is recommended. Etch the enamel surface for 15 – 20 seconds, rinse for 10 seconds and dry in a water and oil-free air stream.

Place LC FLOWFILL SE with the NeedleTip in a thin layer directly onto the prepared areas. With the supplied microbrush work it into the prepared areas under slight pressure for 20s. Remove excess material around the margins with a microbrush.

Cure LC FLOWFILL SE as follows:

Lighter shades (e.g. A1, A2, B1) **20 seconds**

Darker shades (e.g. A3) **30 seconds**

Hold the light emission window as close as possible to the filling material.

Check the occlusion and correct, if necessary.

Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burs.

Polish to high gloss with polishing discs or polishing brushes.

Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

Additional Notes/Warnings

- The curing times are for Halogen curing lights with a light intensity of min. 500 mW/cm² or LED curing lights with a light intensity of min.1000 mW/cm² and a wavelength range 400-500 nm.
- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of LC FLOWFILL SE.
- Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.
- Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.
- Contact with eyes may cause severe eye damage. Wear eye protection. In case of contact with eyes rinse immediately with plenty of water and seek medical advice.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany

Warranty

First Scientific Dental GmbH warrants this product will be free from defects in material and manufacture. First Scientific Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and First Scientific Dental GmbH's sole obligation shall be repair or replacement of the First Scientific Dental GmbH product.

Limitation of Liability

Except where prohibited by law, First Scientific Dental GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!

For dental use only!