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LC IONOLINE

Radiopaque Resin Based Glassionomer Liner

Instructions for Use (EN)

LC IONOLINE is a light cure, radiopaque resin-based glassionomer liner for light cure composites.

Its field of application is:

- lining of cavities for following filling procedures
- lining of cavities to improve marginal behaviour

Side effects

In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calcium hydroxide material.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the material. Neither store the material in proximity of eugenol containing products, nor let the material come into contact with materials containing eugenol.

Application

Clean and dry the cavity.

Screw off the cap of the syringe and substitute it by one of the supplied NeedleTips. Squeeze out the first small drop on a pad to ensure safe handling of LC IONOLINE.

Before applying LC IONOLINE paste use of a radiopaque calcium hydroxide liner, e.g. CALCIUM HYDROXIDE LINER in areas in approximity of the pulp (pulp protection) is recommended. Dentin should be primed with a suitable dentin priming agent (e.g. LC ETCHBOND X1).

Apply then LC IONOLINE in a thin layer until the desired thickness is reached. Light cure LC IONOLINE for **40 seconds**. Do not apply to the cavity edges and remove any paste sticking to them. Then proceed with light cure composite filling.

After use discard the NeedleTip and attach the cap to the syringe again.

NeedleTips are for single use only due to hygienic reasons!

Additional Notes

- Do not use any resin to adjust viscosity of LC IONOLINE.
- Unpolymerized material may have an irritating effect and may lead to a sensitizing reaction against methacrylates.
- Avoid contact with skin, mucous membrane and eyes.
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

Storage

Do not store above 25 °C (77 °F). Do not use after expiry date.

Warranty

First Scientific Dental GmbH warrants this product will be free from defects in material and manufacture. First Scientific Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and First Scientific Dental GmbH's sole obligation shall be repair or replacement of the First Scientific Dental GmbH product.

Limitation of Liability

Except where prohibited by law, First Scientific Dental GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!

For dental use only!