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Endodontic Repair Cement

Instructions for Use (EN)

MTA UNIVERSAL is an endodontic repair cement. MTA UNIVERSAL powder is consisting of very fine hydrophilic particles of several mineral oxides. After contact with MTA UNIVERSAL liquid it forms a gel that hardens to an impermeable barrier.

1. Indications

- Repair of root perforations during root canal therapy
- Root end filling (retrograde)
- Pulp capping
- Root end filling (orthograde)

2. Contraindications

Not known

3. Side effects

Not known

4. Dispensing and Mixing

The **powder/liquid ratio** is **2.6/1.0**. This can be obtained by mixing **1 level (blue) scoop of powder and 2 drops of liquid**.

If a thinner or firmer consistency is desired, the mixing ration can be modified slightly:

Mixing ratio (by weight)	2:1	2,6:1	3:1
Working time (at 23 °C/74 °F)	3:00 min	2:00 min	1:00 min

For root end filling (orthograde) (see 5.4) the recommended mixing ratio is 2:1 (by weight). This can be obtained by mixing **3 level (blue) scoops of powder and 8 drops of liquid**. The mixed material gives a sufficient amount for apexification that can be applied optimally with a suitable application device into the root canal.

For mixing of MTA UNIVERSAL use a mixing pad that is impervious to water or a glass block of suitable dimension.

For accurate dispensing of MTA UNIVERSAL powder shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder with the mixing spatula and carry it onto the mixing pad.

For dispensing of MTA UNIVERSAL liquid turn the bottle vertically with the tip about 5cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop at a time. If any bubbles are present, lightly tap the bottle with the fingers holding it. **Discount under-sized drops** that contain bubbles and are obviously not full-sized. **Discount over-sized drops**, usually resulting from holding the bottle too close to the mixing pad or squeezing the bottle too hard and/or for too long.

Use a small spatula to rapidly mix all the cement powder in portions into the liquid. The mixed cement should be thixotropic and have a homogeneous consistency. Total **mixing time is 30 seconds**.

If desired, a more rigid consistency can be achieved by adding some more powder to the mixture, a more creamy consistency is attained by adding some liquid.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

5. Application

5.1. Repair of root perforations

Place rubber dam and clean the root canal system using intra-canal instruments and irrigate with NaOCl. Dry the root canal with paper points and isolate the perforation.

Fill the apical canal space up to the perforation completely with a suitable root canal filling material.

Mix MTA UNIVERSAL as described under point 4.

Apply MTA UNIVERSAL with suitable instruments into the perforation site and condense it.

Check the position of MTA UNIVERSAL in the root canal by an X-ray. If an adequate barrier has not been created, rinse MTA UNIVERSAL out of the canal and repeat the procedure.

Remove excess moisture with a damp cotton pellet or a paper point.

Place a damp cotton pellet in the access to the root canal and apply a temporary filling material.

Alternatively seal the access preparation with a suitable root canal filling material and seal the cavity with a tight filling.

Both options can be done not before 5 minutes after placement of the MTA UNIVERSAL.

MTA UNIVERSAL repair material remains as a permanent part of the root canal filling.

5.2. Root End Filling (retrograde)

Create an access to the root-end and resect the root with a surgical bur.

Use an ultrasonic tip to prepare a class I root-end cavity preparation to a depth of 3-5 mm.

Isolate the area and dry the root end cavity with paper points. Achieve hemostasis with suitable methods.

Mix MTA UNIVERSAL as described under point 4.

Apply MTA UNIVERSAL with suitable instruments and condense it using a small plugger.

Remove excess cement and clean the surface of the root with a moist piece of gauze. Confirm placement of the MTA-universal repair material with an X-ray. The MTA UNIVERSAL repair material remains as a permanent part of the root canal filling.

5.3. Pulp Capping

Place rubber dam and prepare the cavity outline. If caries is present, remove it. Rinse cavity and exposed pulpal areas with a suitable disinfectant.

Mix MTA UNIVERSAL as described under point 4.

With a suitable instrument apply a small amount of MTA UNIVERSAL over the exposed pulp and remove excess moisture with a dry cotton pellet.

Not before 5 minutes after application of MTA UNIVERSAL place a small amount of a flowable light cure liner (e.g. LC IONOLINE) and light cure.

Apply a suitable bonding agent (e.g. LC ETCHBOND X1) according to the corresponding instructions to the remaining cavity walls.

Place a light cure composite according to the instructions and light cure.

Pulp vitality and status should be checked at regular intervals.

5.4. Root End Filling (orthograde)

Place rubber dam and clean the root canal system using intra-canal instruments and irrigate with NaOCl. Dry the root canal with paper points.

For disinfection place calcium hydroxide paste in the root canal for one week. Seal the access opening with a temporary filling material.

Mix MTA UNIVERSAL as described under point 4..

With a suitable instrument apply a small amount of MTA UNIVERSAL into the apical region and condense it. Create a 3 – 5 mm barrier of MTA UNIVERSAL.

Check the position of MTA UNIVERSAL by an X-ray. If an adequate barrier has not been created, rinse MTA UNIVERSAL out of the canal and repeat the procedure.

Remove excess moisture with a damp cotton pellet or a paper point.

Place a damp cotton pellet in the access to the root canal and apply a temporary filling material.

Alternatively seal the access preparation with a suitable root canal filling material and seal the cavity with a tight filling.

Both options can be done not before 5 minutes after placement of the MTA UNIVERSAL.

MTA UNIVERSAL repair material remains as a permanent part of the root canal filling.

Additional remarks

- In the first hour after application handle the placed MTA cement carefully.
- Intraoral application of MTA UNIVERSAL must be done immediately after mixing to prevent dehydration during setting.
- MTA UNIVERSAL can cause discoloration. Use MTA UNIVERSAL only in the root canal and/or the pulp chamber.

Storage

Store MTA UNIVERSAL at a dry place at 10 – 25 °C (50 °F – 77 °F). **Do not store below 10°C (50 °C)!** Do not use after expiry date.

Warranty

First Scientific Dental GmbH warrants this product will be free from defects in material and manufacture. First Scientific Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and First Scientific Dental GmbH's sole obligation shall be repair or replacement of the First Scientific Dental GmbH product.

Limitation of Liability

Except where prohibited by law, First Scientific Dental GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!

For dental use only!